

Immunofluorescent Assays (IFA) (IVD Medical Devices Class IIb & III – EU Council Directive 98/79/EC)

(Product literature available upon request)

Catalog No.	Description	EDMA No.
IFA- Chagas	Indirect fluorescence assay (IFA) complete kit for diagnosis of Chagas disease(American trypanosomiasis).	15 05 10 02
IFA-Toxoplasmosis	Indirect fluorescence assay (IFA) complete kit for diagnosis of Toxoplasmosis.	15 05 01 00
IFA-Hidatidosis	Indirect fluorescence assay (IFA) complete kit for diagnosis of echinococcosis or hydatid disease.	15 01 03 05
IFA- Syphilis	Indirect fluorescence assay (IFA) complete kit for diagnosis of syphilis.	15 01 03 01
IFA - Triquinosis	Indirect fluorescence assay (IFA) complete kit for diagnosis of trichinosis or trichinellosis or trichinellosis.	15 05 10 01
DFA – Chlamydia T	Direct fluorescence assay (DFA) complete kit for diagnosis of Chlamydia trachomatis antigens in urethral and endocervical swabs samples .	15 01 01 01
IFA – anti DNA	Indirect fluorescence assay (IFA) complete kit for detection of anti native, double stranded DNA using Crith. luciliae.	13 07 01 03
IFA – ANA HEp-2	Indirect fluorescence assay (IFA) complete kit for detection of anti nuclear antibodies using HEp-2 cell	12 10 01 00
IFA-ANA-AMA-ASMA-APCA	Indirect fluorescence assay (IFA) complete kit for detectionof anti nuclear, anti mitochondrial, anti smooth muscle and anti parietal cells antibodies using a	12 10 01 10

	combination of rat liver, kidney and stomach tissues	
IFA – ANA	Indirect fluorescence assay (IFA) complete kit for the detection of anti nuclear , using rat liver tissue	12 10 01 01
IFA – AMA	Indirect fluorescence assay (IFA) complete kit for detection of anti mitochondrial antibodies using rat kidney tissue	12 10 01 01
IFA – ASMA	Indirect fluorescence assay (IFA) complete kit for detection of anti smooth muscle antibodies using rat stomach tissue	12 10 01
IFA – APCA	Indirect fluorescence assay (IFA) complete kit for detection of anti parietal cells antibodies using rat stomach tissue	12 10 01
IFA – AMA-ASMA-APCA	Indirect fluorescence assay (IFA) complete kit for detection of anti mitochondrial, anti smooth muscle and anti parietal cells antibodies using a combination of rat kidney and stomach tissues	12 10 01
IFA – LKM	Indirect fluorescence assay (IFA) complete kit for detection of anti liver kidney microsomes (LKM) antibodies using a combination of rat liver, kidney and stomach tissues.	12 10 01
IFA – EMA IgA	Indirect fluorescence assay (IFA) complete kit For detection of anti endomysial antibodies, IgA class using monkey esophagus	12 10 01