

UROFAST URINE REAGENT TEST STRIPS AND READERS

(Product literature available upon request)

Catalog No.	Description
	UROFAST URINE REAGENT TEST STRIPS
051100	ONE-PARAMETER Blood Glucose Protein
052100	TWO-PARAMETERS Glucose/Ketone Glucose/Protein
053100	THREE-PARAMETERS Glucose/Protein/pH
054100	FOUR-PARAMETERS Glucose/Protein/Ph/Blood Glucose/Protein/pH/Specific Gravity
055100	FIVE-PARAMETERS Blood/Ketone/Glucose/Protein/pH Blood/Urobilinogen/Glucose/Protein/pH
056100	SIX-PARAMETERS Protein/pH/Blood/Ketone/Bilirubin/Glucose
057100	SEVEN-PARAMETERS Urobilinogen/Protein/pH/Blood/Ketone/Bilirubin/Glucose
058100	EIGHT-PARAMETERS Glucose/Protein/pH/Ketone/Urobilinogen/Bilirubin/Blood/Nitrite
059100	NINE-PARAMETERS Glucose/Protein/pH/Ketone/UrobilinogenBilirubin/Blood/Nitrite/Specific Gravity
060100	TEN PARAMETERS Glucose/Protein/pH/Bilirubin/Blood/KetoneUrobilinogen/Ketone/Nitrite/Specific Gravity /Leukocyte
061100	ELEVEN-PARAMETERS Glucose/Protein/pH/Bilirubin/Blood/Ketone/Urobilinogen/Ketone/Nitrite/Specific Gravity/Leukocyte/Ascorbic Acid

UROFAST URINE ANALYZER

UROFAST-100 (ideal for the physician's office or small lab) is a semiautomated and CLIA-waived 60 tests per hour urine strip reader with 1000 test memory. Its space-saving design has patient ID, veterinary test and LIS capability, flip screen, and 110/220 volt compatibility.
