

**BioFast Infectious Disease Rapid Tests** (IVD Medical Devices Class II – EU Council Directive 98/79/EC)

(Product literature available upon request)

Product Name	Description	Specimen	Sensitivity	GMDN CODE
HBsAg	Hepatitis B Virus Surface Antigen	Strip/Cassette Serum/Plasma	1 ng/ml	30723
HBsAb	Hepatitis B Virus Surface Antibody	Strip/Cassette Serum/Plasma	20 mIU/ml	30725
HBeAg	Hepatitis B Virus e Antigen	Strip/Cassette Serum/Plasma	98%	30734
HBeAb	Hepatitis B Virus e Antibody	Strip/Cassette Serum/Plasma	96%	30737
HbCAb	Hepatitis B Virus core Antibody	Strp/Cassette Serum/Plasma	96%	30731
HCV	Hepatitis C Virus Antibody	Strip/Cassette Serum/Plasma	99%	30831
Syph	Syphilis Card Test	Cassette Serum/Plasma		30828
CMV	CMV Card Test	Cassette Serum/Plas		30834
LS	Leptospira	Cassette Serum/Plasma		30716
HIV 1/2	Human Immunodeficiency Virus 1/2 Antibody	Strip/Cassette Serum/Plasma/W. Blood	99%	30769
HIV -Trio (3-line)	Human Immunodeficiency Virus 1/2 Antibody	Strip/Cassette Serum/Plasma/W. Blood	99%	30773
Hep/STD	HBsAg/HCV/HIV Combo	Cassette Serum/Plasma		30778
TP	Treponema pallidum Antibody	Strip/Cassette Serum/Plasma	99%	30828
NG	Neisseria Gonorrhoeae	Strip/Cassette Serum/Plasma	10 <sup>4</sup> cfu/ml	38851
HPAg	Helicobacter Pylori Antigen.	Strip/Cassette Feccal Extract	10 <sup>4</sup> cfu/ml	30825
HPAb	Helicobacter Pylori Antibody	Strip/Cassette Serum/Plasma	95%	30826
M-pf	Malaria Antigen ( pf)	Cassette W.Blood	50/parasite/ml	
M-pf//pv	Malaria Antigen ( pf/pv)	Cassette W.Blood	50 parasite/ml	
M-pf-pan	Malaria Antigen (pf/pan)	Cassette W.Blood	50 parasite/ml	
CV - O1	Cholera Vibrio O1	Strip/Cassette Feccal/Water	10 <sup>4</sup> cfu/ml	
CV-O139	Cholera Vibrio O139	Strip/Cassette Feccal/Water	10 <sup>4</sup> cfu/ml	
CV-Ogawa	Cholera O1 Ogawa	Strip/Cassette Feccal/Water	10 <sup>4</sup> cfu/ml	

CTA	Chlamydia T Antigen	Cassette	Vag Swab		30676	
LT	Lyme Test	Cassette	Blood			
LS	Leptosira Card Test	Cassette	Feccal/Water/W Blood		30716	
YP- F1	Plague F1 Antigen Y.pestis	Strip/Cassette	Feccal Extract	10 <sup>4</sup> cfu/ml		
TB	Tuberculosis (TB)	Cassette	Serum/Plasma	95%		
DENG-Ab	Dengue Virus -IgG/IgM Antibody	Strip/Cassette	Serum/Plasma/W.Blood		39519	
DENG- Ag	Dengue Virus Antigen Combo	Cassette	Serum/W Blood		39519	
DENG NS1	Dengue NS1 Antigen	Strip/Cassette	Serum/W Blood	1 ng/ml	39519	
ADE	Adenovirus	Strip/Cassette	Feccal Extract	10 particle/ml	30832	
ROT	Rota Virus	Strp/Cassette	FecalExtract	10 particle/ml	30832	
ROT-ADE	Rotavirus+ Adenovirus ( 3 lines)	Strp/Cassette	Feccal Extract	10 particle/ml		
PROCAL	Procalcitonin Test	Strp/Cassette	Serum/Plasma			
STREP A STREP B	Streptococcus Group A Streptococcus Group B	Cassette	Throat Swab	10 <sup>5</sup> cfu/ml	30826/30827	
CHLAM-T	Chlamydia Trachomatis	Cassette	Vag/Ure Swab	10 <sup>5</sup> cfu	30678	
S-T-Ag	Salmonella Typhi Antigen	Strp/Antigen	Serum/PlasmaWBloodFecca		30707	
S-T-PT-Ag	Salmonella Typhi/Para Typhi Antigen	Strip/Cassette	Serum/PlasmaWBlood/Fec		30707	
CHK	Chikungunya IgG/IgM	Strip/Cassette	Serum/Plasma			
CHAG	Chagas Antibody	Strp/Cassette	Serum/Plasma		30820	
Z	Zika Virus	Strip/Cassette	Serum/Plasma			

## Other Rapid Diagnostics Tests (IVD Medical Devices Class II – EU Council Directive 98/79/EC)

(Product literature available upon request)

Product Name	Description	Specimen	GMDN Code
<b>TUMOR MARKERS</b>			
CEA	Carcinoembryonic Antigen	Cassette Serum/Plasma	17293
AFP	Alpha Fetoprotein	Cassette Serum/Plasma	3029
PSA	Prostate Specific Antigen Total	Cassette Serum/Plasma	30511
HPL	HPL Card Test	Cassette Serum/Plasma	30514
FOB-g	Fecal Occult Blood –Gastric fluid	Cassette Fecal Extract	17410
FFOB	Fecal Occult Blood and Fecal Transferrin Combo test	Cassette Fecal Extract	
FT	Fecal Transferrin	Cassette Serum/Plas	17402
<b>CARDIAC MARKERS</b>			
CK-MB	Creatine Kinase Isoenzymes	Cassette Serum/W Blood	30261
MYO	Myoglobin	Cassette Serum/W Blood	30508
TROP	Troponin I	Cassette Serum/W Blood	30509
CK/MY/TR	CK-MB/Myoglobin/Troponin I Combo	Cassette Serum	
FABP	Fatty Acid Binding Protein	Cassette Serum/W Blood	30262
CRP	C- Reactive Protein (Qualitative & Semi Quantitative)	Cassette Serum/W Blood	17290
<b>FERTILITY</b>			
HCG	Foremost Pregnancy Test Urine	Strip/Cassette Urine	30513
HCG- S/U	Foremost Pregnancy Test Combo	Strip/Cassette Urine/Serum	30513
HCG-MS	Foremost Midstream Pregnancy Test	Cassette Urine	30513
HCG- Early	Foremost Early Detection Pregnancy Test	Cassette Urine	30513
LH	Foremost Luteinizing Hormone Ovulation Test	Strip/Cassette Urine	30515
FSH	Foremost Follicle Stimulating Hormone	Strip/Cassette Urine	30512
MFT	Foremost Male Fertility Test	Cassette Sperm	
MFT-OCT	Foremost Male Fertility Test- OCT	Cassette Sperm	
MFT-POCT	Foremost Male Fertility Sperm Test- POCT	Cassette Sperm	